

In addition to the apparatus and fittings found in a chemistry laboratory, each candidate will require the following:

1. about 150 cm³ of solution A labelled solution A ✓
2. about 150 cm³ of solution B labelled solution B ✓
3. about 80 cm³ of solution C labelled solution C ✓
4. one pipette 25.0 ml
5. one pipette filler
6. one volumetric flask (250.0 ml)
7. four labels
8. about 500 cm³ of distilled water
9. one burette 50.0 ml
10. three conical flasks
11. one 10 ml measuring cylinder
12. one 100ml measuring cylinder
13. two boiling tubes
14. one thermometer -10° C to 110° C
15. about 0.5 g of solid E supplied in a stoppered container
16. six clean dry test-tubes
17. about 0.1 g of solid F supplied in a stoppered container
18. about 0.5 g of solid G supplied in a stoppered container
19. pH chart 1-14; and universal indicator solution supplied with a dropper
20. one 100ml beaker
21. one **metallic** spatula
22. one clean dropper

Access to

1. Phenolphthalein indicator supplied with a dropper ✓
2. 2 M sulphuric (VI) acid supplied with a dropper ✓
3. 2 M sodium hydroxide supplied with a dropper ✓
4. 0.5 M potassium iodide supplied with a dropper ✓
5. Bromine water supplied with a dropper ✓
6. Acidified potassium manganate (VII) supplied with a dropper ✓
7. Bunsen burner. ✓

PREPARATIONS

1. Solution A is prepared by taking 180.0 cm³ of concentrated hydrochloric acid (specific gravity 1.18) adding it to 600 cm³ of distilled water in a 1 litre volumetric flask and diluting it to the mark. Label this solution as **solution A**.
2. Solution B is prepared by dissolving 80.0 g of sodium hydroxide pellets in 800 cm³ of distilled water and diluting it to the mark. Label it as **solution B**.
3. Solution C is prepared by dissolving 25g of solid C in 600 cm³ of distilled water and diluting it to the mark. Label this as **solution C**.
4. Bromine water is prepared by taking 1 cm³ of liquid bromine and dissolving it in 100 cm³ of distilled water in a fume cupboard. This must be freshly prepared and supplied in a dropper bottle.
5. Acidified potassium manganate (VII) is prepared by dissolving 3.16g of solid potassium manganate (VII) in about 600 cm³ of 2 M sulphuric (VI) acid and adding distilled water to make 1 litre.

N.B: Solids C, E, F and G will be supplied by the KNEC.

1st Session

Candidates' Index Numbers

.....

.....

.....

.....

.....

.....

.....

Question 1: Draw and complete table I as per the question paper.

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

Question 1: Draw and complete table II as per the question paper.

.....

.....

.....

.....

.....

This form **MUST** be completed and returned in the envelope containing the scripts.

**Kenya Certificate of Secondary Education
Oct./Nov. 2010**

REPORT ON CHEMISTRY PRACTICAL 233/3

The supervisor is required to:

- (a) give the names and index numbers of candidates in each session.
- (b) ask the teacher in charge of chemistry to perform the experiments in question 1 **(NOT IN THE SAME ROOM WITH CANDIDATES)** and write the results in the tables as per the question paper.
- (c) give the details of any difficulties experienced by particular candidates writing their names and index numbers on page 8.

Details required include:

- (i) difficulties caused by faulty apparatus
- (ii) accidents caused by chemicals or apparatus
- (iii) physical handicaps e.g colour blindness supported by medical evidence.

Name of Supervisor:

Centre Number:

Signature:

Date:

Name of the teacher who performed the experiments:

TSC No:

Signature: **Date:**

Seating plan.

www.atikaschool.com

2nd Session

Candidate's Index Numbers

.....

.....

.....

.....

.....

.....

Question 1: Draw and complete table I as per the question paper.

.....

.....

.....

.....

.....

.....

.....

.....

Question 1: Draw and complete table II as per the question paper.

.....

.....

.....

.....

.....

.....

Seating plan.